

PT. BINA BARUTAMA SEJAHTERA

Specific Gravity Of Semi Solid Bituminous Materials

BB 501 ASTM D.70 / AASHTO T.228







BB 501.2



Specific Gravity Of Semi **Solid Bituminous Materials**

Туре	Parts Name	Specification	Qty	
BB 501.1	Picnometer	Hubbard type 24 ml capacity	1 P	С
BB 501.2	Picnometer	Hubbard type 25 ml capacity	1 P	С
BG 034/4	Thermometer	Range 0 – 200 C	1 P	С

Supplied without:

BG.011/1 Hot Plate

BG.045/3 Asbestos Wire Gauze

BG.027/1 Glass Funnel Weight / volume relationships are useful for shipping purposes and for the calculation of mix proportions during mix design. the test method determines the ratio of a given volume of binder to the weight of an equal volume of water at a given temperature. Usually 25 C Hubbard Carmick type pycnometer are used, being particularly suitable as they combine ease of filling and cleaning with long term stability.